



86079

## Neissers Stain A (Methylene Blue)

Part A

### Specifications

Appearance (Colour)	Dark blue
Appearance (Form)	Liquid
Suitability as a stain for microscopy	Passes test

### Other Information

#### Applications

Solution for diphtheria staining in Microbiological, suitable for diagnostic assay manufacturing.  
Suitable for hematology, histology.

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	24 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
125 ml	38229090 (GST 12%)

### Available Packages

125 ml

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com