



86133

Furfural extrapure AR, ACS, 99%

Part B

CAS : 98-01-1

Molecular Formula : C₅H₄O₂

Molecular Weight : 96.09

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Pale yellow
Appearance (Form)	Liquid
Assay (GC)	min. 99%
Density (g/ml) @ 20°C	1.158-1.162
Refractive Index (20°C)	1.521-1.525
Boiling Range	161-163°C
Non Volatile Matter	max. 0.5%
Acidity (C ₅ H ₄ O ₃)	max. 0.5%
Water (KF)	max. 0.05%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1199 , IMCO Class No.:6.1(3) , Packing Group:II
HSN Code	
500 ml	29122990 (GST 18%)

Available Packages

500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com