



86263

N-Methyl-2-Pyrrolidone (NMP) extrapure AR, 99.5%

Part A

CAS : 872-50-4

Molecular Formula : C₅H₉NO

Molecular Weight : 99.13

Specifications

| | |
|-------------------------|-------------|
| Appearance (Clarity) | Clear |
| Appearance (Colour) | Colourless |
| Appearance (Form) | Liquid |
| Assay (GC) | min. 99.5% |
| Density (g/ml) @ 20°C | 1.029-1.031 |
| Refractive Index (20°C) | 1.470-1.471 |
| Boiling Range | 201-202°C |
| Non Volatile Matter | max. 0.001% |
| Water (KF) | max. 0.1% |

General Information

| | |
|---------------------|---------------------------------------|
| Storage | 25 to 40°C (Room Temperature) |
| Shelf Life | 60 Months |
| IMDG Identification | Not Regulated for Transport (Non-Haz) |
| HSN Code | |
| 2500 ml | 29142990 (GST 18%) |
| 500 ml | 29142990 (GST 18%) |
| Type of Packing | |
| 2500 ml | Glass Bottle |
| 500 ml | Glass Bottle |

Available Packages

500 ml

2500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com