



86543

## Jenners Stain (Eosin-methylene blue)

Part A

CAS : 62851-42-7

### Specifications

Appearance (Colour)	Greenish black
Appearance (Form)	Powder
Solubility (Turbidity) 0.1% solution in Methanol	Clear
Solubility (Colour) 0.1% solution in Methanol	Blue solution with green fluorescence
Loss on drying	max. 10%
A 1% methanol	
@524nm	min. 600
@654nm	min. 1100
Suitability as a stain for microscopy	Passes test

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	32049000 (GST 18%)
25 Gms	32049000 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle
25 Gms	Glass Bottle

### Available Packages

25 Gms  
100 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com