



86570

Triethylenetetramine (mix. of isomers) pure, 97%

Part A

CAS : 112-24-3

Molecular Formula : C₆H₁₈N₄

Molecular Weight : 146.24

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless to pale yellow
Appearance (Form)	Liquid
Assay (GC) (mixture of Isomers)	min. 97%
Density (g/ml) @ 20°C	0.970-0.980
Refractive Index (20°C)	1.498-1.499

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:2259 , IMCO Class No.:8 , Packing Group:II
HSN Code	
2500 ml	29212990 (GST 18%)
500 ml	29212990 (GST 18%)
Type of Packing	
500 ml	Glass Bottle

Available Packages

500 ml

2500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com