



86633

Silver Iodate extrapure, 99%

Part A

CAS : 7783-97-3

Molecular Formula : AgIO₃

Molecular Weight : 282.77

Specifications

Appearance (Colour)	White to off-white
Appearance (Form)	Powder
Assay	min. 99%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1479 , IMCO Class No.:5.1 , Packing Group:II
HSN Code	
10 Gms	28299030 (GST 18%)
100 Gms	28299030 (GST 18%)
25 Gms	28299030 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle

Available Packages

10 Gms
25 Gms
100 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com