



86701

Lithium Aluminium Hydride extrapure, 95%

Part B

CAS : 16853-85-3

Molecular Formula : LiAlH₄

Molecular Weight : 37.95

Specifications

| | |
|---------------------|---------------|
| Appearance (Colour) | White to grey |
| Appearance (Form) | Powder |
| Assay | min. 95% |

General Information

| | |
|---------------------|---|
| Storage | 25 to 40°C (Room Temperature) |
| Shelf Life | 60 Months |
| IMDG Identification | UN No.:1410 , IMCO Class No.:4.3 , Packing Group:I |
| HSN Code | |
| 100 Gms | 28500010 (GST 18%) |
| Type of Packing | |
| 100 Gms | Glass Bottle |

Available Packages

100 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com