



86715

Paraffin Wax Cubes (Type 1 - 56-58) extrapure

Part A

CAS : 8002-74-2

Specifications

| | |
|---------------------|------------------------------|
| Appearance (Colour) | White to Off white |
| Appearance (Form) | Cubes (~ 20 -40mm) or Chunks |
| Melting Point | 56 - 58°C |

General Information

| | |
|---------------------|---------------------------------------|
| Storage | 25 to 40°C (Room Temperature) |
| Shelf Life | 60 Months |
| IMDG Identification | Not Regulated for Transport (Non-Haz) |
| HSN Code | |
| 1000 Gms | 27129090 (GST 18%) |
| 500 Gms | 27129090 (GST 18%) |
| Type of Packing | |
| 1000 Gms | Plastic Jar |
| 500 Gms | Plastic Bottle |

Available Packages

500 Gms
1 Kg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com