



86936

L-Amino Acids Kit

Part B

Specifications

Assay	min. 99% (All L-Amino acids in the kit)
TLC test	Single spot passes

Other Information

Notes

All 22 L-Amino acids - 1g each, in the kit, have assay min.99% and TLC test shows single spot. Amino acids contained are,

Amino acids contained are,
L-Amino acids plus Glycine, L-Alanine,
L-Arginine Hydrochloride, L-Asparagine,
L-Aspartic acid, L-Cysteine free base,
L-Cystine,
L-Glutamic acid, L-Glutamine, Glycine,
L-Histidine Hydrochloride, 4-Hydroxy-
L-Proline,
L-Isoleucine, L-Leucine, L-Lysine
Hydrochloride,
L-Methionine, L-Phenylalanine, L-
Proline,
L-Serine, L-Threonine, L-Tryptophan,
L-Tyrosine, L-Valine.

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Kit	38229090 (GST 12%)

Available Packages

1 Kit

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com