



87326

10X Tris-Glycine pH-8.3 Tank Buffer for molecular biology

Part B

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
pH	8.3
Sterility (Bacterial count)	Nil
Suitability for Electrophoresis	Passes test
DNase, RNase, Protease	Not detected

General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	36 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1000 ml	29420090 (GST 18%)
200 ml	29420090 (GST 18%)

Available Packages

200 ml

1000 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com