



87367

## Cupric Nitrate Trihydrate extrapure, 98%

Part A

CAS : 10031-43-3

Molecular Formula :  $\text{Cu}(\text{NO}_3)_2 \cdot 3\text{H}_2\text{O}$

Molecular Weight : 241.60

### Specifications

Appearance (Colour)	Blue
Appearance (Form)	Hygroscopic compound
Assay	min. 98%
Chloride (Cl)	max. 0.005%
Sulphate (SO <sub>4</sub> )	max. 0.02%
Iron (Fe)	max. 0.02%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1477 , IMCO Class No.:5.1 , Packing Group:II
HSN Code	
250 Gms	28342990 (GST 18%)
500 Gms	28342990 (GST 18%)
Type of Packing	
250 Gms	Glass Bottle
500 Gms	Glass Bottle

### Available Packages

250 Gms  
500 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com