



87767

Buffer capsules pH 4 (10 capsules)

Part L

Specifications

pH

4.0 +/- 0.05
(Dissolve contents of 1 capsule
& dilute upto 100 ml distilled water
to prepare solution of pH 4.0 at 20°C)

Other Information

Notes

Dissolve contents of 1 capsule & dilute
upto 100 ml distilled water to prepare
solution of pH 4.0 at 20°C

General Information

Storage

25 to 40°C (Room Temperature)

Shelf Life

60 Months

IMDG Identification

Not Regulated for Transport (Non-Haz)

HSN Code

10 caps

38229090 (GST 12%)

Available Packages

10 caps

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com