



88001

1,1,3,3-Tetramethyldisiloxane (TMDS) extrapure, 98%

Part B

CAS : 3277-26-7

Molecular Formula : C₄H₁₄OSi₂

Molecular Weight : 134.32

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Assay (GC)	min. 98%
Density (g/ml) @ 20°C	0.757-0.758
Refractive Index (20°C)	1.370-1.371
Boiling Range	70-71°C

General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
HSN Code	
100 ml	29319090 (GST 18%)
25 ml	29319090 (GST 18%)

Available Packages

25 ml
100 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com