



88147

## DPX Mountant for histology

Part A

### Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Viscous liquid
Density (g/ml) @ 20°C	0.920 - 0.921
Refractive Index (20°C)	~ 1.520

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1307 , IMCO Class No.:3 , Packing Group:III
HSN Code	
1000 ml	29420090 (GST 18%)
250 ml	29420090 (GST 18%)
500 ml	29420090 (GST 18%)

### Available Packages

250 ml

500 ml

1000 ml

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com