



88383

Tetrabutylammonium Hydrogen Sulphate (TBAHS) for HPLC, 99.5%

Part A

CAS : 32503-27-8

Molecular Formula : C₁₆H₃₇NO₄S

Molecular Weight : 339.54

Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 10% aq. solution	Clear
Solubility (Colour) 10% aq. solution	Colourless
Assay	min. 99.5%
Melting Point	171 - 174°C
Sulphated Ash	max. 0.05%
Loss on drying	max. 0.5%
Absorbance (A) of 0.005 M aq. solution in a 1cm cell v/s H ₂ O	
@200nm	max. 0.15
@220nm	max. 0.05
@250 nm	max. 0.01

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	29239000 (GST 18%)
25 Gms	29239000 (GST 18%)
500 Gms	29239000 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle
25 Gms	Glass Bottle

Available Packages

25 Gms

100 Gms

500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com
