



88459

## FMOC-L-Arginine extrapure, 98%

Part B

CAS : 91000-69-0

Molecular Formula : C<sub>21</sub>H<sub>24</sub>N<sub>4</sub>O<sub>4</sub>

Molecular Weight : 396.45

### Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Assay	min. 98%

### General Information

Storage	0 to 4°C (Freeze)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Gms	29379020 (GST 18%)
100 Gms	29379020 (GST 18%)
25 Gms	29379020 (GST 18%)
5 Gms	29379020 (GST 18%)

### Available Packages

1 Gms

---

5 Gms

---

25 Gms

---

100 Gms

---

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com