



88639

## 2,4-Dinitrophenylhydrazine extrapure AR, 99%

Part B

CAS : 119-26-6

Molecular Formula : C<sub>6</sub>H<sub>6</sub>N<sub>4</sub>O<sub>4</sub>

Molecular Weight : 198.14

### Specifications

Appearance (Colour)	Orange red
Appearance (Form)	Crystalline compound moistened with water
Assay (On dry basis)	min. 99%
Sulphated Ash	max. 0.02%
Loss on drying	~30%
Insoluble matter in sulphuric acid	max. 0.01%

### Other Information

Notes	(Forbidden by Air)
-------	--------------------

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:3380 , IMCO Class No.:4.1 , Packing Group:I
HSN Code	
100 Gms	29280090 (GST 18%)
25 Gms	29280090 (GST 18%)
500 Gms	29280090 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle
25 Gms	Glass Bottle
500 Gms	Glass Bottle

### Available Packages

25 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srlchem.com, website www.srlchem.com