



88726

L-Threonine for cell culture, 99%, Endotoxin (BET) 1EU/mg

Part B1

CAS : 72-19-5

Molecular Formula : C₄H₉NO₃

Molecular Weight : 119.12

Specifications

Appearance (Colour)	White to off white
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 10% solution 2N HCl	Clear
Solubility (Colour) 10% solution 2N HCl	Colourless to pale yellow
Assay (NT)	min. 99%
Specific rotation [?]	-29° to -27° (C=2,water)
Loss on drying	max. 0.5%
Sulphated Ash	max. 0.1%
Chloride (Cl)	max. 0.01%
Sulphate (SO ₄)	max. 0.01%
Iron (Fe)	max. 0.001%
Heavy Metals (Pb)	max. 0.001%
Ammonium (NH ₄)	max. 0.01%
Homogeneity by TLC	Passes test
Suitability for Cell Culture	Passes test
Bacterial Endotoxin Test (BET) (EU/mg = IU/mg)	1 EU/mg

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	29224990 (GST 18%)
25 Gms	29224990 (GST 18%)
Type of Packing	
100 Gms	Plastic Bottle
25 Gms	Plastic Bottle

Available Packages

25 Gms

100 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com
