



88917

## Tungsten Nanopowder, 50nm, 99.9%

Part C

CAS : 7440-33-7

Atomic Formula : W

Atomic Weight : 183.84

### Specifications

Assay	min. 99.9%
APS	50nm

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
HSN Code	
1 Gms	81011000 (GST 18%)
5 Gms	81011000 (GST 18%)

### Available Packages

1 Gms  
5 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com