



89049

Bismuth (III) Oxide extrapure, 99%

Part A

CAS : 1304-76-3

Molecular Formula : Bi<sub>2</sub>O<sub>3</sub>

Molecular Weight : 465.96

Specifications

Appearance (Colour)	Yellow
Appearance (Form)	Powder
Assay	min. 99%
Nitrate (NO <sub>3</sub> )	max. 0.03%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	28259090 (GST 18%)
500 Gms	28259090 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle
500 Gms	Glass Bottle

Available Packages

100 Gms
500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com