



89141

Lithium Lactate (specially purified) extrapure, 99%

Part B

CAS : 867-55-0

Molecular Formula : C₃H₅O₃Li

Molecular Weight : 96.01

Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline compound
Solubility (Turbidity) 10% aq. solution	Clear
Solubility (Colour) 10% aq. solution	Colourless
Assay (NT)	min. 99%
Loss on drying	max. 0.5%
Chloride (Cl)	max. 0.005%
Sulphate (SO ₄)	max. 0.01%
Iron (Fe)	max. 0.005%
Heavy Metals (Pb)	max. 0.003%
Arsenic (As)	max. 0.0001%
Acidity (Lactic acid)	max. 0.1%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	29181190 (GST 18%)
500 Gms	29181190 (GST 18%)
Type of Packing	
100 Gms	Plastic Bottle

Available Packages

100 Gms

500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com