



89159

Xylene extrapure AR, ACS, ExiPlus, Multi-Compendial, 99.5%

Part A

CAS : 1330-20-7

Molecular Formula : C8H10

Molecular Weight : 106.17

Specifications

ExiPlus	meets compendial specs of USP
Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Colour (APHA)	max. 10
Total Assay (GC, C8H10 isomers + Ethylbenzene)	min. 99.5%
Ethyl benzene (GC)	max. 25%
Density (g/ml) @ 20°C	0.860-0.865
Refractive Index (20°C)	1.497-1.499
Boiling Range	136-140°C
Non Volatile Matter	max. 0.001%
Substance Discoloured by H2SO4	Passes test
Acidity (HCl)	max. 0.001%
Iron (Fe)	max. 0.00001%
Heavy Metals (Pb)	max. 0.00001%
Sulphur Compound (CS2)	max. 0.0003%
Thiophene homologues (C4H4S)	max. 0.0001%
Benzene (GC)	max. 0.1%
Toluene (GC)	max. 0.5%
Water (KF)	max. 0.05%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1307 , IMCO Class No.:3 , Packing Group:III
HSN Code	
2500 ml	29024400 (GST 18%)
500 ml	29024400 (GST 18%)
Type of Packing	

2500 ml

Glass Bottle

Available Packages

500 ml

2500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srlchem.com, website www.srlchem.com