



89184

D-Arginine (free base) extrapure CHR, 99%

Part B

CAS : 157-06-2

Molecular Formula : C₆H₁₄N₄O₂

Molecular Weight : 174.20

Specifications

| | |
|---------------------|--------------------|
| Appearance (Colour) | White |
| Appearance (Form) | Crystalline powder |
| Assay (NT) | min. 99% |
| Loss on drying | max. 0.5% |
| Homogeneity by TLC | Passes test |

General Information

| | |
|---------------------|---------------------------------------|
| Storage | 25 to 40°C (Room Temperature) |
| Shelf Life | 60 Months |
| IMDG Identification | Not Regulated for Transport (Non-Haz) |
| HSN Code | |
| 1 Gms | 29224990 (GST 18%) |
| 5 Gms | 29224990 (GST 18%) |
| Type of Packing | |
| 25 Gms | Glass Bottle |

Available Packages

1 Gms

5 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com