



89299

Diisopropylether Dried, 99%, water 0.0025%

Part A

CAS : 108-20-3

Molecular Formula : C6H14O

Molecular Weight : 102.18

Specifications		
Appearance (Clarity)	Clear	
Appearance (Colour)	Colourless	
Appearance (Form)	Liquid	
Assay (GC)	min. 99%	
Density (g/ml) @ 20°C	0.722-0.724	
Refractive Index (20°C)	1.368-1.369	
Boiling Range	67-69°C	
Acidity (CH3COOH)	max .0.005%	
Peroxide (H2O2)	max. 0.005%	
Water (KF)	max. 0.0025%	

Other Information		
Notes	Packed with Dr. Dried Seal which are designed with corrosion free Crown Cap & with specially developed coloured septa. This enable repetitive needle puncturing easy and securing the consistency of the product from outside atmosphere.	

General Information		
Storage	25 to 40°C (Room Temperature)	
Shelf Life	60 Months	
IMDG Identification	UN No.:1159 , IMCO Class No.:3 , Packing Group:II	
HSN Code		
1000 ml	29091990 (GST 18%)	
500 ml	29091990 (GST 18%)	
Type of Packing		
1000 ml	Glass Bottle	

Available Packages

500 ml

1000 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: [info@srlchem.com](mailto:info@srlchem.com), website [www.srlchem.com](http://www.srlchem.com)