



89320

Soya Lecithin (30%) BactoBio for bacteriology

Part D

CAS : 8002-43-5

Specifications

Appearance (Colour)	Cream to yellow
Appearance (Form)	Free flowing powder
pH (25°C)	6.5 - 7.5 (2%)
Acid value	= 36.0%
Moisture & Volatile	Maximum 1.5%
Phosphatidylcholine (PC) & Lyso Phosphatidylcholine (LPC)	Minimum 30% (HPLC Method)
Total Aerobic microbial count	Maximum 1000 cfu/gm
Total Yeast and Mould count	Maximum 100 cfu/gm
Salmonella/25gm	Absent
Escherichia coli	Absent
Coliforms	10 cfu/gm

Other Information

Applications

A natural phospholipid made from non-GMO soya beans.

General Information

Storage	8 to 25°C (Cool & Dry Area)
Shelf Life	36 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	29232010 (GST 18%)
500 Gms	29232010 (GST 18%)
Type of Packing	
500 Gms	Plastic Bottle

Available Packages

100 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com