



89429

## Brij-35® (Main Component) (Polyoxyethylene Lauryl Ether, Brij-L23)

Part A

CAS : 9002-92-0

### Specifications

Appearance (Colour)	White to pale yellow
Appearance (Form)	Waxy solid
Solubility (Turbidity) 10% aq. solution	Clear to slight hazy
Acid value	max. 3.0
Water (KF)	max. 3.0%

### Other Information

Notes	A non-ionic detergent for use in various protein methods.
-------	---

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
250 Gms	34024200 (GST 18%)
500 Gms	34024200 (GST 18%)
Type of Packing	
250 Gms	Plastic Bottle
500 Gms	Plastic Bottle

### Available Packages

250 Gms  
500 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com