



89635

Cyclohexanone extrapure AR, 99.5%

Part A

CAS : 108-94-1

Molecular Formula : C6H10O

Molecular Weight : 98.15

Specifications		
Appearance (Clarity)	Clear	
Appearance (Colour)	Colourless	
Appearance (Form)	Liquid	
Colour (APHA)	max. 10	
Assay (GC)	min. 99.5%	
Density (g/ml) @ 20°C	0.945-0.947	
Refractive Index (20°C)	1.450-1.451	
Boiling Range	154-156°C	
Non Volatile Matter	max. 0.02%	
Acidity (CH3COOH)	max. 0.02%	
Iron (Fe)	max. 0.00005%	
Heavy Metals (Pb)	max. 0.00005%	
Cyclohexane (GC)	max. 0.001%	
Cyclohexanol (GC)	max. 0.2%	
Copper (Cu)	max. 0.00001%	
Water (KF)	max. 0.05%	

General Information		
Storage	25 to 40°C (Room Temperature)	
Shelf Life	60 Months	
IMDG Identification	UN No.:1915 , IMCO Class No.:3 , Packing Group:III	
HSN Code		
2500 ml	29142200 (GST 18%)	
500 ml	29142200 (GST 18%)	
Type of Packing		
500 ml	Glass Bottle	

Available Packages

500 ml

2500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srlchem.com, website www.srlchem.com