



90251

Graphene Oxide Film (Super Paper) (GOFSP),20000nm, 98%

Part C

CAS : 7782-42-5

Molecular Formula : C

Molecular Weight : 12.01

Specifications

Appearance (Colour)	Black
Appearance (Form)	Circular paper film
Assay	min. 98% C
Thickness	20 μ
Diameter	40 mm

Other Information

Application

Applications
Preparation of Graphene, Solar energy, Graphene semiconductor chips, Conductive graphene film, Graphene computer memory, Biomaterials, Transparent conductive coatings, etc.

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 piece	25049090 (GST 5%)

Available Packages

1 piece

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com