



90365

Phenylphosphate Disodium Salt Dihydrate (High Purity) extrapure AR, 98%

Part B

CAS : 66778-08-3

Molecular Formula : C₆H₅Na₂O₄P·2H₂O

Molecular Weight : 254.09

Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 10% aq. solution	Clear
Solubility (Colour) 10% aq. solution	Colourless
Assay	min. 98%
pH (10% aq. solution)	9.0 - 11.0
Chloride (Cl)	max. 0.01%
Iron (Fe)	max. 0.0005%
Heavy Metals (Pb)	max. 0.0002%
Phosphate (P)	max. 0.1%
Free Phenol	max. 0.002%
Water (KF)	~ 15%

Other Information

Application	For detection of phosphatase in milk
-------------	--------------------------------------

General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	29199090 (GST 18%)
25 Gms	29199090 (GST 18%)
250 Gms	29199090 (GST 18%)
Type of Packing	
100 Gms	Plastic Bottle
25 Gms	Plastic Bottle
250 Gms	Plastic Bottle

Available Packages

25 Gms

100 Gms

250 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com
