



90464

## Pectinase ex. Aspergillus Niger, 3U/mg powder

Part B

CAS : 9032-75-1

### Specifications

Appearance (Colour)	Off-white to pale brown
Appearance (Form)	Powder
Activity	min. 3.5 units/mg protein
Unit Definition	One unit of pectinase activity causes an increase of 0.01 absorbance per minute with Pectin as substrate at pH 5.8 at 30 °C in a reaction volume of 3.10 mL.
Enzyme solution for activity determination	Prepared in citrate buffer, pH 5.8.

### General Information

Storage	-20 °C (Blue/Dry Ice)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1000 Units	35079050 (GST 18%)

### Available Packages

1000 Units

---

5000 Units

---

10000 Units

---

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com