



90544

Trichloroacetic Acid extrapure, 99%

Part A

CAS : 76-03-9

Molecular Formula : C₂HO₂Cl₃

Molecular Weight : 163.39

Specifications

Appearance (Colour)	White
Appearance (Form)	Very hygroscopic crystals
Solubility (Turbidity) 100% aq. solution	Clear
Solubility (Colour) 100% aq. solution	Colourless
Assay (by Neutralization)	min. 99%
Melting Point	55 - 58°C
Sulphated Ash	max. 0.02%
Chloride (Cl)	max. 0.02%
Sulphate (SO ₄)	max. 0.02%
Iron (Fe)	max. 0.002%

General Information

Storage	8 to 25°C (Cool & Dry Area)
Shelf Life	60 Months
IMDG Identification	UN No.:1839 , IMCO Class No.:8 , Packing Group:II
HSN Code	
100 Gms	29154030 (GST 18%)
500 Gms	29154030 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle
500 Gms	Glass Bottle

Available Packages

100 Gms
500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com