



90691

Neutral Red extrapure AR, 60%

Part A

CAS : 553-24-2

Molecular Formula : C₁₅H₁₇CIN₄

Molecular Weight : 288.78

Specifications

Appearance (Colour)	Black to brown to dark green
Appearance (Form)	Powder
Solubility (Turbidity) 0.1% solution in 50% ethanol + 2 ml acetic acid?	Clear
Loss on drying	max. 12%
Dye content	min. 60 %
Absorbance (A) of 1% soln in 50% Ethanol + 2ml acetic acid in a 1cm cell @540nm	min. 900

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
10 Gms	32049000 (GST 18%)
100 Gms	32049000 (GST 18%)
25 Gms	32049000 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle
25 Gms	Glass Bottle

Available Packages

10 Gms
25 Gms
100 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com