



90820

## p-Anisic Acid extrapure, 98%

Part B

CAS : 100-09-4

Molecular Formula : C<sub>8</sub>H<sub>8</sub>O<sub>3</sub>

Molecular Weight : 152.15

### Specifications

Appearance (Colour)	White to off white to pale cream
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 5% hot solution in Methanol	Clear
Assay (NaOH titration)	min. 98%
Melting Point	180 - 186°C
Sulphated Ash	max. 0.2%
Water (KF)	max. 0.5%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	29163190 (GST 18%)
500 Gms	29163190 (GST 18%)
Type of Packing	
100 Gms	Plastic Bottle
500 Gms	Plastic Bottle

### Available Packages

100 Gms

500 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com