



91186

Xylan from Beechwood

Part B

CAS : 9014-63-5

Specifications

Appearance (Colour)	Pale beige to Brown
Appearance (Form)	Powder
Solubility (Turbidity) 1% solution in 5% NaOH	Clear to slightly opalescent
Solubility (Colour) 1% solution in 5% NaOH	Brown
Specific rotation [?] @25°C	-40° to -80° (C=0.5%, 5%NaOH) (on dry basis)
Loss on drying	max. 5%
Suitability as a substrate for enzyme Xylanase	Passes

General Information

Storage	8 to 25°C (Cool & Dry Area)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
10 Gms	29400000 (GST 18%)
25 Gms	29400000 (GST 18%)
5 Gms	29400000 (GST 18%)

Available Packages

5 Gms

10 Gms

25 Gms

100 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com