



91229

Zinc Metal Granules, 99%, 3-8 mm

Part A

CAS : 7440-66-6

Molecular Formula : Zn

Molecular Weight : 65.39

Specifications

Appearance (Colour)	Silver or Grey
Appearance (Form)	Granules or needles
Assay (Complexometric)	min 99.0%
Mesh size (3 to 8 mm)	min.80%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
HSN Code	
500 Gms	79031000 (GST 18%)

Available Packages

500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com