



91482

## Titanium Ethoxide (ET) extrapure, 33-35% TiO<sub>2</sub>

Part A

CAS : 3087-36-3

Molecular Formula : C<sub>8</sub>H<sub>20</sub>O<sub>4</sub>Ti

Molecular Weight : 228.11

### Specifications

Appearance (Colour)	Colourless to pale yellow
Appearance (Form)	Liquid
Assay (TiO <sub>2</sub> content)	33 - 35%
Density (g/ml) @ 20°C	1.088 - 1.089

### Other Information

Notes

For Nanotechnology, Organic Synthesis, Catalysis, Polymerization of Epoxy's, Thin Films and Cross-Linking.

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1993 , IMCO Class No.:3 , Packing Group:III
HSN Code	
100 ml	29051490 (GST 18%)
250 ml	29051490 (GST 18%)
Type of Packing	
250 ml	Glass Bottle

### Available Packages

100 ml

250 ml

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com