



91496

Tetramethylammonium Hydroxide 10% aq. solution pure

Part A

CAS : 75-59-2

Molecular Formula : C₄H₁₃NO

Molecular Weight : 91.15

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless to pale yellow tinted
Appearance (Form)	Liquid
Assay (T)	~10%
Halides (Cl)	max. 0.02%
Sulphated Ash	max. 0.1%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1835 , IMCO Class No.:8 , Packing Group:III
HSN Code	
250 ml	38229090 (GST 12%)
50 ml	38229090 (GST 12%)
Type of Packing	
250 ml	Plastic Bottle
50 ml	Plastic Bottle

Available Packages

50 ml
250 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com