



91695

4-Benzyloxy-Phenylboronic Acid extrapure, 97%

Part B

CAS : 146631-00-7

Molecular Formula : C13H13BO3

Molecular Weight : 228.05

Specifications

Appearance (Colour)	Pale yellow
Appearance (Form)	Crystalline powder
Assay (NaOH titration)	min. 97%
Melting Point	198 - 200°C

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
HSN Code	
1 Gms	29319090 (GST 18%)
25 Gms	29319090 (GST 18%)
5 Gms	29319090 (GST 18%)

Available Packages

1 Gms
5 Gms
25 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srlchem.com, website www.srlchem.com