



91740

## Lithium Iodide Monohydrate pure, 98%

Part B

CAS : 17023-24-4

Molecular Formula : LiI.H<sub>2</sub>O

Molecular Weight : 151.86

### Specifications

Appearance (Colour)	White to yellow or brown
Appearance (Form)	Crystalline compound / flakes
Assay (On dry basis)	min. 98%
Sulphate (SO <sub>4</sub> )	max. 0.03%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	28276090 (GST 18%)
25 Gms	28276090 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle

### Available Packages

25 Gms  
100 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com