



91837

## N,N-Diisopropylethylamine (DIPEA) extrapure, 98%

Part A

CAS : 7087-68-5

Molecular Formula : C<sub>8</sub>H<sub>19</sub>N

Molecular Weight : 129.25

### Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless to pale yellow
Appearance (Form)	Liquid
Assay (GC)	min. 98%
Density (g/ml) @ 20°C	0.741-0.743
Refractive Index (20°C)	1.413-1.414
Boiling Range	126-128°C
Water (KF)	max. 0.3%

### Other Information

Notes	Suitable for peptide synthesis. (Forbidden by Air)
-------	---

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:3384 , IMCO Class No.:6.1(3) , Packing Group:I (Forbidden by Air)
HSN Code	
100 ml	29211990 (GST 18%)
1000 ml	29211990 (GST 18%)
250 ml	29211990 (GST 18%)
Type of Packing	
100 ml	Glass Bottle
1000 ml	Glass Bottle
250 ml	Glass Bottle

## Available Packages

100 ml

250 ml

1000 ml

---

### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: [info@srchem.com](mailto:info@srchem.com), website [www.srchem.com](http://www.srchem.com)

---