



92035

## 40% Acrylamide / Bis-acrylamide Mix Solution (Ratio: 37.5:1)

Part A

### General Information

Storage	2 to 8°C (Refrigerate)
IMDG Identification	UN No.:3426 , IMCO Class No.:6.1 , Packing Group:III
HSN Code	
100 ml	38229090 (GST 12%)
500 ml	38229090 (GST 12%)

### Available Packages

100 ml

500 ml

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com