



92101

Tannic Acid extrapure AR, ACS, ExiPlus, Multi-Compendial

Part A

CAS : 1401-55-4

Molecular Formula : C₇₆H₅₂O₄₆

Molecular Weight : 1701.22

Specifications

ExiPlus	meets compendial specs of BP
Appearance (Colour)	Off-white to yellowish brown
Appearance (Form)	Powder
Identification Test	Passes test
Loss on drying	max.12%
Residue after ignition	max. 0.5%
Zinc (Zn)	max. 0.005 %
Heavy Metals (Pb)	max. 0.003 %
Sugars, Dextrin	Passes test

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	32019030 (GST 18%)
250 Gms	32019030 (GST 18%)
500 Gms	32019030 (GST 18%)
Type of Packing	
100 Gms	Plastic Bottle

Available Packages

100 Gms

250 Gms

500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com