



92196

## Biphenyl extrapure, 99%

Part B

CAS : 92-52-4

Molecular Formula : C<sub>12</sub>H<sub>10</sub>

Molecular Weight : 154.21

### Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline powder or crystals or flakes
Solubility (Turbidity) 5% solution in Methanol	Clear
Solubility (Colour) 5% solution in methanol	Colourless
Assay (GC)	min. 99%
Melting Point	69 - 72°C
Sulphated Ash	max. 0.05%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:3077 , IMCO Class No.:9 , Packing Group:III
HSN Code	
100 Gms	29029090 (GST 18%)
1 Kg	29029090 (GST 18%)
Type of Packing	
1 Kg	Plastic Jar
100 Gms	Glass Bottle

### Available Packages

100 Gms

1 Kg

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com