



92391

## Hexamine extrapure AR, 99.5%

Part A

CAS : 100-97-0

Molecular Formula : C<sub>6</sub>H<sub>12</sub>N<sub>4</sub>

Molecular Weight : 140.19

### Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 10% aq. solution	Clear
Solubility (Colour) 10% aq. solution	Colourless
Assay (NaOH Titration, On dry basis)	min. 99.5 %
pH (10% aq. solution)	8.5 - 9.5
Chloride (Cl)	max. 0.002%
Sulphate (SO <sub>4</sub> )	max. 0.005%
Iron (Fe)	max. 0.0005%
Sulphated Ash	max. 0.05%
Heavy Metals (Pb)	max. 0.0001%
Copper (Cu)	max. 0.0001%
Ammonium (NH <sub>4</sub> )	max. 0.005%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1328 , IMCO Class No.:4.1 , Packing Group:III
HSN Code	
5 Kg	29212910 (GST 18%)
500 Gms	29212910 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle
500 Gms	Glass Bottle

### Available Packages

500 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com