



92665

Pyridine AR, ACS, ExiPlus, Multi-Compendial, 99.5%

Part A

CAS : 110-86-1

Molecular Formula : C₅H₅N

Molecular Weight : 79.10

Specifications

ExiPlus	meets compendial specs of USP, BP, Ph.Eur
Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Solubility with water	No turbidity
Assay (GC)	min. 99.5%
Density (g/ml) @ 20°C	0.980-0.983
Refractive Index (20°C)	1.508-1.510
Boiling Range	114-116°C
Chloride (Cl)	max. 0.0005%
Sulphate (SO ₄)	max. 0.0005%
Ammonia (NH ₃)	max. 0.002%
Copper (Cu)	max. 0.0005%
Reducing Substances	Passes test
Non Volatile Matter	max. 0.002%
Water (KF)	max. 0.1%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1282 , IMCO Class No.:3 , Packing Group:II
HSN Code	
250 ml	29333100 (GST 18%)
2500 ml	29333100 (GST 18%)
500 ml	29333100 (GST 18%)
Type of Packing	
500 ml	Glass Bottle

Available Packages

250 ml

500 ml

2500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com
