



92822

p-Nitrophenyl-β-D-Glucuronide (PNPG (Gus)) extrapure, 99%

Part B

CAS : 10344-94-2

Specifications

Appearance (Colour)	White to off white to pale yellow
Appearance (Form)	Crystalline powder
Assay (HPLC)	min. 99.0%

General Information

Storage	0 to 4°C (Freeze)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Mg	29095090 (GST 18%)
50 Mg	29095090 (GST 18%)
500 Mg	29095090 (GST 18%)
Type of Packing	
100 Mg	Glass Bottle
50 Mg	Glass Bottle

Available Packages

50 Mg
100 Mg
500 Mg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com