



92875

Biochemical Test Kit - I (Pre-Weighed Media)

Part D

Other Information

Applications

Used for biochemical identification of gram negative enteric bacteria based on tests like Carbohydrate utilization, H₂S production, Phenylalanine deaminase and Urease production.

Composition

(55010)Peptone Water	15.00 gm/lt.
(55720)Triple Sugar Iron Agar	64.63 gm/lt.
(92335)Phenylalanine Agar	26.00 gm/lt.
(18339)Urea Agar Base	24.00 gm/lt.

General Information

Storage	8 to 25°C (Cool & Dry Area)
Shelf Life	36 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Kit	38229090 (GST 12%)

Available Packages

1 Kit

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com