



92901

FITC-Dextran 70000 ex. Leuconostoc Sp.

Part B

CAS : 60842-46-8

Specifications

Appearance (Colour)	Yellow to orange to yellowish orange
Appearance (Form)	Powder or crystals
Average Molecular weight	63 - 77 KDa
Loss on drying	max. 10%
Degree of substitution	0.001 - 0.008

General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Mg	32129090 (GST 18%)
250 Mg	32129090 (GST 18%)
Type of Packing	
100 Mg	Glass Bottle
250 Mg	Glass Bottle

Available Packages

100 Mg
250 Mg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com