



92996

Barium Chloride Dihydrate ACS, Exiplus, Multi-Compendial, 99%

Part A

CAS : 10326-27-9

Molecular Formula : BaCl₂.2H₂O

Molecular Weight : 244.28

Specifications

ExiPlus	meets compendial specs of USP, BP, Ph.Eur
Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Solubility (Colour) 10% aq. solution	Colourless
Solubility (Turbidity) 10% aq. solution	Clear
Insoluble matter	max. 0.005%
Assay (ex Cl)	min. 99%
pH (5% aq. solution)	5.0 - 8.0
Iron (Fe)	max. 0.0001%
Loss on drying	14.0 - 16.0 %
Copper (Cu)	max. 0.0002%
Heavy Metals (Pb)	max. 0.0002%
Nitrogen Compounds (N)	max. 0.002%
Potassium (K)	max. 0.0025%
Sodium (Na)	max. 0.005%
Calcium (Ca)	max. 0.005%
Strontium (Sr)	max. 0.1%
Oxidising substances	max. 0.005%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1564 , IMCO Class No.:6.1 , Packing Group:III
HSN Code	
5 Kg	28273990 (GST 18%)
500 Gms	28273990 (GST 18%)
Type of Packing	
500 Gms	Plastic Bottle

Available Packages

500 Gms

5 Kg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com
