



93119

Mercuric Thiocyanate extrapure AR, ExiPlus, Multi-Compendial, 98%

Part A

CAS : 592-85-8

Molecular Formula : Hg(SCN)₂

Molecular Weight : 316.75

Specifications

ExiPlus	meets compendial specs of Ph.Eur
Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Assay (Complexometric)	min. 98%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1646 , IMCO Class No.:6.1 , Packing Group:II
HSN Code	
100 Gms	28521000 (GST 18%)
25 Gms	28521000 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle
25 Gms	Glass Bottle

Available Packages

25 Gms
100 Gms
250 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srlchem.com, website www.srlchem.com